

# Update

Dear Practice Nurses and Doctors,

We have changed our collection requirements for Zinc, Copper and other trace elements.

## As of 27<sup>th</sup> July 2023:

For Zinc and Copper testing we are now utilising  
a new Trace Elements tube  
(Dark Blue with Red stripe on label)  
which contains a clot activator coating.



For all other trace elements,  
we are utilising EDTA tube (Lavender).



## Doctors Order Request

Surgery :	
Address	
Contact Name:	Phone:
Date:	Received: (office use)

### Doctor's request for stores **GOODS TO BE USED FOR 4CYTE PATHOLOGY COLLECTIONS ONLY.**

Under Medicare Australia guidelines, collection equipment is only supplied if used for pathology collection. Fortnightly ordering is recommended for efficient stock rotation. Please allow 5 working days for delivery.

DESCRIPTION (UNIT OF ISSUE)	QTY	MAX	DESCRIPTION (UNIT OF ISSUE)	QTY	MAX
<b>VACCUUM BLOOD COLLECTION</b>			<b>URINE/STERILE COLLECTION</b>		
Barrel for multisample needles (Bag of 250)		1	Specimen Jar (Yellow top) (Each)		50
Needle: 21G X 1 1/2" multisample (Box 100)		1	24hr Plain Urine Collection Container: Nil preservative (Each)		2
Needle: 22G x 1 1/2" multisample (Box 100)		1	24hr Acid Urine Collection Container (Each)		2
Needle: 21G Butterfly Green (Each)		10	24hr Urine Collection Funnel (Each)		2
Needle: 23G Butterfly Blue (Each)		10	Paediatric Urine Bag (Each)		10
Tube: EDTA (Lavender top) 4ml (Each)		100	<b>GYNAECOLOGICAL AND CYTOLOGICAL COLLECTION</b>		
Tube: SST (Gold top) 8.5ml (Each)		100	ThinPrep® (Tray 25)		2
Tube: Sodium Citrate (Light Blue top) 2.7ml (Each)		100	Rovers/Wallach Cervex Sampler (Pack 25)		2
Tube: EDTA (Blue top) Trace Element 6ml (Each)		6	Endo-Cervical Brush (Pack 10)		2
Tube: Fluoride Oxalate (Grey top) 2.0ml (Each)		100	FLOQ Swab Self Collect CST (Each)		10
Tube: Lithium Heparin (Green top) 4ml (Each)		6	Speculum – EOS Small (Each)		20
Tube: EDTA Gel PPT (White/Pearl top) 5ml (Each)		6	Speculum – EOS Medium (Each)		20
Tube: Plain Clot (Red top) 6ml (Each)		6	Speculum – EOS Large (Each)		20
Microtainer Tube: EDTA (Lavender) 1ml (Each)		1	Speculum – VIKING Small (Each)		20
Microtainer Tube: SST (Gold) 1ml (Each)		1	Speculum – VIKING Medium (Each)		20
Microtainer Tube: Fluoride Oxalate (Grey) 1ml (Each)		1	Speculum – VIKING Large (Each)		20
Tube: Citrate (Light Blue) 1.8ml (Each)		1	Cotton Applicator: Single tip: 15cm (Pack 100)		1
Blood Culture Bottle – Aerobic (Blue Cap) (Each)		1	<b>HISTOPATHOLOGY</b>		
Blood Culture Bottle – Anaerobic (Purple Cap) (Each)		1	Histology Container: 25ml with formalin (Tray 24)		1
Blood Culture Bottle – Paeds (Pink Cap) (Each)		1	Histology Container: 70ml with formalin (Tray 10)		1
Kit: Percept NIPT (VCGS)		4	Biopsiblade® (Box 25)		1
Kit: Quantiferon TB Gold Plus (Set of 4 tubes)		5	Biopsy Punch: 3mm with plunger (Box 20)		1
<b>NON-EVACUATED BLOOD COLLECTION</b>			Biopsy Punch: 4mm with plunger (Box 20)		1
Syringe: 5ml without needle (Box 100)		1	Biopsy Punch: 6mm (Box 20)		1
Syringe: 10ml without needle (Box 100)		1	Biopsy Punch: 8mm (Box 20)		1
Needle: 21G x 1 1/2" (Box 100)		1	<b>STATIONERY</b>		
Needle: 23G x 1 1/4" (Box 100)		1	A4 Pathology Request Form with labels (Pack 100)		5
<b>FAECES COLLECTION</b>			Pathology Request Pad (Carbon) (Book)		2
Faeces Jar (Brown top) (Each)		20	A4 Maternal Serum Screen Form (Pad) VCGS		1
Kit: FOB (Set of 3 Jars)		10	A4 Percept NIPT Form (Pad) VCGS		1
<b>SWAB COLLECTION</b>			A4 Prepair Genetic Carrier Screening Form (Pad) VCGS		1
eSwab Regular Collection (Blue cap) (Each)		20	Doctor's Request for Stores (Sheet) (Each)		10
eSwab Minitip (NPG/Urethral) (Orange Cap) (Each)		20	Swab and Brush Guide		6
<b>GENERAL</b>			Tube Guide		6
Clear Specimen Bag (Pack 100)		5	<b>MISCELLANEOUS</b>		
Red Specimen Bag (Pack 100)		2			
Glucose Load Drink 75gm (Each)		5			
Alcohol Swab (Pack 100)		3			
Urea Breath Test Kits C14 (Each)		5			